

## Dräger SSK 5000

The Draeger SSK 5000 kit is used for sampling for drugs on or from surfaces. When used in combination with the Dräger DrugTest® 5000, the sample is transferred to a test cassette and screened for possible drug contamination. For further testing, the collected sample can be stored and transported.



The need to monitor for substance abuse or misuse means that, more and more frequently law enforcement officers are required to check for potential surface contamination and/or unknown substances on-site. This needs to be carried out quickly and reliably. The routine drug analysis performed by forensic laboratories, which are equipped for this task, is not always possible in terms of time or logistics.

With the new Dräger SSK 5000, used in conjunction with the Dräger DrugTest® 5000, it is now possible to reliably test drug contaminated surfaces or drug-suspicious substances quickly and on-site.

The Dräger SSK 5000 consists of a tube of buffer solution and, in the cap, an integral swab for sampling. The swab is used to wipe surfaces and, subsequently, test for drugs and their metabolites via the Dräger DrugTest® 5000 System.

Once transferred into the buffer solution using the Dräger SSK 5000, the surface sample or unknown substance can be stored and transported safely in the sample container. A special screw cap enables the sample container to be converted into a dropping bottle. The sample solution can then be separated or reduced for subsequent analysis using droplet dosing techniques.




The Dräger SSK 5000: for easy, fast, and safe surface sampling, storage and transportation.


### ORDER INFORMATION

Item	Order No.
Dräger DrugTest® 5000 Analyzer	83 19 900
Dräger DrugTest® 5000 Test Kit 6-drug panel (20 per box)	83 19 965
<b>Accessories</b>	
Dräger SSK 5000 Device for surface samples (20 per box)	83 19 910


**DRÄGER SSK 5000 - SAMPLING PROCEDURE:**



Wipe approx. 15.5 in<sup>2</sup> of the surface with the moisturized swab.




Insert the swab into the SSK 5000 buffer tube, seal and shake for approx. 30 seconds.




Using the Test Kit cassette: Open the tube cassette, upside down into the tube. Turn the cassette right side up and wait until the sample indicator has turned blue.


or



Alternate option for 3a- Using the dropping bottle: Slowly apply 12 to 13 drops onto the test cartridge collector. Important: do not use more than 13 drops.



Seal the container and, if necessary, store for further inspection.



Perform an analysis with the Draeger DrugTest® 5000 Analyzer according to user manual instructions.

**Storage and Stability**

The Draeger SSK 5000 must be stored at a temperature between 39 °F to 86 °F.

[www.draeger.com](http://www.draeger.com)

**USA**

Draeger Safety Diagnostics, Inc.  
4040 West Royal Lane, Suite 136  
Irving, TX 75063

**USA**

Tel +1 972 929-1100  
Fax +1 972 929-1105

[info@dsgi.biz](mailto:info@dsgi.biz)

**CANADA**

Draeger Safety Canada Ltd.  
7555 Danbro Crescent  
Mississauga, ON L5N 6P9  
Canada

Tel +1 905 821-8988  
Fax +1 905 821-2565

[sales.canada@draeger.com](mailto:sales.canada@draeger.com)

**MEXICO**

Draeger Safety S.A. de C.V.  
Av. Pe-uelas No. 5, Bodega No. 37  
Fraccionamiento Industrial San Pedrito  
Queretaro, Qro C.P. 76148

Mexico  
Tel +1 52 442 246 1113  
Fax +1 52 442 246 1114

[ventas.mexico@draeger.com](mailto:ventas.mexico@draeger.com)

**GERMANY**

Dräger Safety AG & Co. KGaA  
Revalstrasse 1  
23560 Lübeck

**Germany**

Tel +49 451 882-0  
Fax +49 451 882-2080